

## Settings

### Display

---

..: You can modify the interface language: *english* or *french*. Just click on your choice.

**Note**-This change requires the restart of *Alimenthèque™*.

..: You can customize the software by inserting your logo. The selected logo will appear on all the printed reports.

**Note**-The software creates a link towards the picture file desired. If the picture file is moved, the link will be lost and will simply need to be recreated.

..: You can determine the default standard that will appear in every new recipe created, as well as the default standard that will appear in every new dietary evaluation.

You can also select the proteins standard's unit (*g/day* or *g/kg/day*).

..: You can also determine units of some nutrients: energy (kcal ou kJ), vitamin D ( $\mu\text{g}$  ou IU), vitamin E (mg ou IU) et niacin (mg ou NE).